

## CRYOLAB 25 inc accessories

For dispensing into small  
dewars and flasks.

The following is included:

- Cryolab 25.  
(25 litres of LN<sub>2</sub>)
- Withdrawal Device inc  
Swan neck  
& Phase separator.
- Roller base.

When ordering request  
Part number **9901130**

1. Easy handling of the 25 litre Dewar weighing 26 kg
2. Below PED, so no yearly safety checks required- saving money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

The following part numbers are included in this system:  
9901062 - 9701056 - 9901063



## CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.  
(25 litres of LN<sub>2</sub>)
- Withdrawal Device inc  
a 1.2 Meter hose  
& Phase separator.
- Roller base.  
Ease of movement.

When ordering request  
Part number 9901131

1. Easy handling of the 25 litre Dewar weighing 26kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.



## CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.  
(25 litres of LN<sub>2</sub>)
- Withdrawal Device inc  
a 1.2 Meter hose  
& Phase separator.
- Transport trolley.  
Ease of movement.

When ordering request  
Part number **9901132**



The pouring Trolley is  
1033mm high x 472mm wide.

1. Easy handling of the 25 litre Dewar weighing 26kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.

## CRYOLAB 25 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 25.  
(25 litres of LN2)
- Withdrawal Device  
including Swan neck  
& Phase separator.
- Transport trolley.  
Ease of movement.

When ordering request  
Part number **9901133**

1. Easy handling of the 25 litre Dewar weighing 26 kg
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

The following part numbers are included in this system:  
9901062-9701056-9901064



The pouring Trolley is  
1033mm high x 472mm wide.

## CRYOLAB 35 inc accessories

For dispensing into small  
dewars and flasks.

The following is included:

- Cryolab 35.  
(35 litres of LN<sub>2</sub>)
- Withdrawal Device  
including Swan neck  
& Phase separator.
- Roller base

When ordering request  
Part number **9901134**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with phase separator.

The following part numbers are included in this system:  
9901080 - 9701057A - 9902082



## CRYOLAB 35 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 35.  
(35 litres of LN<sub>2</sub>)
- Withdrawal Device  
including a  
1.2 Meter hose  
& Phase separator.
- Roller base

When ordering request  
Part number **9901135**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. Can be filled by any gas company.
5. Dispense up to 2 litres/minute with limited head pressure.
6. Dispenser fitted with hose and phase separator.

The following part numbers are included in this system:  
9901080 - 9701152A - 9902082



## CRYOLAB 50 inc accessories

For dispensing into small  
dewars and flasks.


The following is included:

- Cryolab 50.  
(50 litres of LN2)
- Withdrawal Device  
including Swan neck  
& Phase separator.
- Roller base

When ordering request  
Part number **9901136**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. 50% the cost of a stainless steel vessel and will keep LN2 longer.
5. Can be filled by any gas company.
6. Dispense up to 2 litres/minute with limited head pressure.
7. Dispenser fitted with phase separator.

The following part numbers are included in this system:  
9901090-9701057-9902082



## CRYOLAB 50 inc accessories

For refrigerator top up.

The following is included:

- Cryolab 50.  
(50 litres of LN2)
- Withdrawal Device  
including a  
1.2 Meter hose  
& Phase separator.
- Roller base

When ordering request  
Part number **9901137**

1. Easy to manoeuvre trolley.
2. Below PED, so no yearly safety checks required- save on money.
3. Aluminium Dewar so a lighter construction than stainless steel vessel and easier to move.
4. 50% the cost of a stainless steel vessel and will keep LN2 longer.
5. Can be filled by any gas company.
6. Dispense up to 2 litres/minute with limited head pressure.
7. Dispenser fitted with hose and phase separator.

The following part numbers are included in this system:  
9901090-9701152-9902082

