

# PRESSURE EQUIPMENT DIRECTIVE 97/23/EC

## THE PRESSURE EQUIPMENT REGULATIONS 1999 (S.I. 2001)

### DECLARATION OF CONFORMITY

We hereby certify that the pressure equipment stipulated below complies with all the relevant provisions of the Pressure Equipment Directive and the National Laws and Regulations adopting this Directive.

Manufacturer	Name: Statebourne Cryogenics Ltd., Address: 18 Parsons Road, Parsons Ind.Est., Washington, Tyne & Wear, England. Post Code: NE37 1EZ.
--------------	---

Vessel	Part Number	PER Category	Approx. Overall Dim.
Cryostor 30	9911025	I	786mm high x 540mm dia
Cryostor 60	9911037	I	900mm high x 620mm dia
Cryostor 75	9911044	II	1051mm high x 620mm dia
Cryostor 90	9911045	II	1803mm high x 620mm dia
Cryostor 120	9911065	II	1122mm high x 720mm dia
Cryostor 180	9911075	II	1378mm high x 720mm dia
Cryostor 200	9911150	II	1445mm high x 750mm x 750mm
Cryostor 240	9911087	II	1645mm high x 690mm x 690mm
Cryostor 300	9911089	II	1439mm high x 818mm dia
Cryostor 300	9911155	III	1550mm high x 900mm x 900mm
Cryostor 400	9911290	III	1660mm high x 820mm x 820mm
Cryostor 500	9911300	III	1913mm high x 820mm x 820mm
Cryostor 500	9911303	III	1937mm high x 820mm x 1120mm

EC Type Certificate No .....  
EC Design Certificate No .....  
EC Certificate of Conformity .....  
Conformity Procedure used. - Modules H1

Notes .....

Harmonized Standards Applied (including parts/clauses) EN 10028-7: 2007  
EN 287-1: 2004  
EN 15614-1:2004


Standards/Specifications Applied (including parts/clauses) - PD5500:2009  
AD Merkblatt A2: 2004  
CGA Pamphlet S1.1 (Section 5.9.1)  
EN 10088-2:1995


Other Community Directives Applied .....

Notified Body	Name: BSI (Notified Body No 0086) Address: Maylands Avenue Hemel Hempstead, England. Post Code: HP2 4SQ.
---------------	---

  
Approved For Engineering

  
Approved For Quality Assurance

  
Approved For Production

  
Managing Director

Issue: 7  
Date:- 08/01/2010