

**STATEBOURNE**

**cryogenics**

# *OD Series*

*Maxi & Mini*  
Industrial Open Dewars

1 Year Vacuum Warranty

Laboratory Open Dewars

Roller Bases

Wide Neck



**5** STATEBOURNE  
cryogenics  
years

*Cryogenic innovation & technology  
you can trust...*

## Industrial Open Dewars Maxi



OD Maxi 80

A Mini and Maxi range of stainless steel Open Dewars for laboratory process, bench top and industrial freezing applications provides capacities from 0.5 litres to 80 litres.

Of exceptionally strong stainless steel construction, these dewars are ideal for use as cold traps, for immersion or simply for temporary liquid storage.



OD Maxi 25



OD Maxi 15



OD Insulating lids

All Maxi models are supplied with insulating lids and are fitted with carry handles. Roller base available for Maxi 15, 25, 80

## Product Features

- OD Series 8 models, from 0.6 litres to 66 litres capacity Open Dewars
- Roller bases available
- Low Cost
- High thermal efficiency
- OPEN DEWAR 1 YEAR VACUUM WARRANTY



Perforated baskets available for the Open Dewar series see Data.

## Laboratory Open Dewars Mini

A Mini range of stainless steel Open Dewars for laboratory process, bench top and industrial freezing applications provides capacities from 0.5 litres to 5.5 litres.

Of exceptionally strong stainless steel construction, these dewars are ideal for use as cold traps, for immersion or simply for temporary liquid storage.

All models (except the OD 0.5) are supplied with insulating lids and swivel carry handles. (The 0.5 Mini is not supplied with a lid or handle.)

OD Mini 3



OD Mini 6



OD Mini 1



OD Mini 0.5

(The 0.5 Mini is not supplied with a lid or handle.)

OD Mini 2

# OD Mini Series

## Technical Data

MODEL	OD Mini 0.5	OD Mini 1	OD Mini 2	OD Mini 3	OD Mini 6
<b>PART NUMBER</b>	<b>9915010</b>	<b>9915020</b>	<b>9915030</b>	<b>9915040</b>	<b>9915050</b>
Volume Litres	0.6	0.9	1.8	2.6	5.5
Int Diameter mm	65	85	100	185	185
External Diameter mm	87	107	122	200	200
Internal Height mm	180	206	285	160	270
External Height mm	204	232	312	190	300
Wt Kg Empty	0.5	0.9	1.3	1.8	2.4
Wt Kg Full	1.0	1.6	2.8	3.9	6.8
Roller Base Part Number	-	-	-	-	-
Perforated Basket	-	-	9915036	-	9915051

# OD Maxi Series

## Technical Data

MODEL	OD Maxi 15	OD Maxi 25	OD Maxi 80
<b>PART NUMBER</b>	<b>9915080</b>	<b>9915090</b>	<b>9915115</b>
Volume Litres	18	31	66
Internal Diameter mm	260	327	406
External Diameter mm	333	385	460
Internal Height mm	370	400	513
External Height mm	565	595	777
Wt Kg Empty	13	15.0	28.0
Wt Kg Full	27.5	40.0	81.0
Roller Base Part Number	9915082	9915092	9915125
Perforated Basket	9915086	9915096	9915117

*Statebourne Cryogenics*  
 18 Parsons Road,  
 Parsons Industrial Estate,  
 Washington,  
 Tyne and Wear  
 NE37 1EZ  
 England.

**Web:** [www.statebourne.com](http://www.statebourne.com)

**Tel:** +44 (0) 191 416 4104

**Fax:** +44 (0) 191 415 0369

**Email:** [info@statebourne.com](mailto:info@statebourne.com)

